

Cell Cytotoxicity Fluorometric Assay Kit

Cat. No. Kit-0196 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Cell Cytotoxicity Assay Kit (Fluorometric) is a non-radioactive and no-wash method assay suitable for high throughput screening of cell proliferation or cytotoxicity against a variety of compounds.
Description	The measurement of mitochondrial dehydrogenases (e.g. LDH) activity is a well-accepted assay to quantify cell numbers and monitor cell viability. Our Cell Cytotoxicity Assay Kit (Fluorometric) provides a fast, simple, accurate and homogeneous assay for the colorimetric or fluorimetric detection of viable cells. This assay is based on the observation that oxidized non-fluorescent blue resazurin is reduced to a red fluorescent dye (resorufin) by accepting an electron from mitochondrial respiratory chain in live cells. The amount of resorufin produced is directly proportional to the number of living cells.
Applications	This Cell Cytotoxicity Assay Kit (Fluorometric) is more sensitive for cell proliferation and cytotoxicity than other assays such as MTT.
Storage	Keep in freezer and protect from light.
Kit Components	Component A: Assay Solution 20 mL
Compatible Sample Types	Cells
Features & Benefits	Non-Radioactive: No special requirements for waste treatment. Continuous: Easily adapted to automation without mixing or separation. Convenient: Formulated to have

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

minimal hands-on time. Wide Applications: Cell proliferation and cytotoxicity. Sensitive
And Accurate: As low as 100 cells can be accurately quantified. Enhanced Value:
Less expensive than the sum of individual components.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA