

Cell Viability NIR Fluorescence Assay Kit

Cat. No. Kit-0202 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Cell Viability Assay Kit (NIR Fluorescence) is a fluorescent assay used to label cells for the research of cellular functions.
Description	The far red fluorescent dye used in the kit is a hydrophobic compound that can easily permeate intact live cells and gets enhanced fluorescence upon entering into live cells. The hydrolysis of the weakly fluorescent substrate by intracellular esterases generates a strongly red fluorescent hydrophilic product that is well-retained in the cell cytoplasm. The red fluorophore generated by the weakly fluorescent substrate has the spectral properties similar to those of Cy5 or Alexa Fluor 647 at Ex/Em = ~640/660 nm (Cy5 or Alexa Fluor 647 are the trademarks of GE Healthcare and Invitrogen respectively). When well excited with the He-Ne Laser at 633 nm, the red fluorophore emits intense red fluorescence at ~660 nm.
Applications	The Cell Viability Assay Kit (NIR Fluorescence) is a set of tool for monitoring cell viability and cellular functions.
Usage	For research use only (RUO)
Storage	Keep in freezer and protect from light.
Kit Components	Component A: CytoCalcein NIR (lyophilized) 2 vials Component B: DMSO 100 µL Component C: Assay Buffer 10 mL x 2
Detection method	Fluorometric

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

**Compatible Sample
Types**

Cells

Features & Benefits

Robust: Higher maximum signal with lower variation across the plate. Convenient:
Formulated to have minimal hands-on time. Rapid Dye Loading: Dye loading at RT for
30 min to 1 hr. Versatile applications: Compatible with many cell lines and targets.

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA