

Fixable Dead Cell Staining Kit (Green Fluorescence)

Cat. No. Kit-0274 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Fixable Dead Cell Staining Kit (Green Fluorescence) is a set of tool used to evaluate the viability of mammalian cells by flow cytometry.
Description	The Fixable Dead Cell Staining Kit (Green Fluorescence) is a set of tools used to evaluate the viability of mammalian cells by flow cytometry. The fluorescent dyes provided in the kits are retained in cells by reacting with cellular components. For viable cells, only the cell-surface amines are available to react with the dyes while for the necrotic cells or the other cells with compromised membranes, the reactive dyes react with cell surface amines and intracellular amines, resulting in more intense fluorescent staining. The difference in fluorescence intensity between the live and dead cell populations is ~100-500 folds and can be completely preserved after fixation.
Usage	For research use only (RUO)
Storage	Keep in freezerProtect from moisture and lightAvoid repeated freeze-thaw cycles
Kit Components	Component A: Reactive Dye 1 vialComponent B: DMSO 200 µl
Compatible Sample Types	Cells

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA