

DNA Damage Quantification Kit

Cat. No. Kit-0296 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	DNA Damage Quantification Kit is used for determining the DNA lesion and repair.
Description	Apurinic/aprimidinic (AP) sites are one of the major types of DNA lesions formed during the course of base excision and repair of oxidized, deaminated or alkylated bases. It has been estimated that about 2×10^5 base lesions are generated per cell per day. The level of AP sites in cells can be a good indicator of DNA lesion and repair against chemical damage and cell aging. The DNA Damage Quantification Kit utilizes the ARP (Aldehyde Reactive Probe) reagent that reacts specifically with an aldehyde group which is the open ring form of the AP sites. After treating DNA containing AP sites with ARP reagents, AP sites are tagged with biotin residues, which can be quantified using avidin-biotin assay followed by a colorimetric detection. The kit provides the necessary reagents for convenient determination of abasic sites in purified DNA sample in 96-well plate format.
Target Species	Mammals
Usage	For research use only (RUO)
Storage	Store the kit at 4°C. Do not freeze.
Kit Components	ARP Solution (10 mM, in DMSO) 0.125 ml TE Buffer 30 ml Glycogen Solution (10 g/l) 0.1 ml 0 ARP-DNA Standard (0.5 g/ml) 0.6 ml 40 ARP-DNA Standard (0.5 g/ml) 0.6 ml DNA Binding Solution 10 ml HRP-Streptavidin 0.15 ml 10X Wash Buffer 30 ml HRP Developer 10 ml 96-well Microplate (8 x 12 strips) 1 plate

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Detection method	Colorimetric
------------------	--------------