

DNMT3B Assay Kit

Cat. No. Kit-0310 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	DNMT3B Assay Kit is use for measuring the amount of DNMT3B.
Description	Epigenetic inactivation of genes plays a critical role in many important human diseases, especially in cancer. A core mechanism for epigenetic inactivation of the genes is methylation of CpG islands in genome DNA. Methylation of CpG islands involves the course in which DNA methyltransferases (Dnmts) transfer a methyl group from S-adenosyl-L-methionine to the fifth carbon position of the cytosines. At least three families of Dnmts have been identified in mammals now: Dnmt1, Dnmt2, and Dnmt3. The Dnmt3 family comprises of three different proteins: Dnmt3A, Dnmt3B, and Dnmt3L. Dnmt3A and Dnmt3B have been demonstrated to methylate both unmethylated and hemimethylated DNA equally and supposed to mediate de novo methylation together with Dnmt1. Increased activation or amount of Dnmt3 is believed to be involved in carcinogenesis, and other genetic and epigenetic diseases. Several methods such as Western blot are
Applications	For measuring Dnmt amount quantitatively from fresh tissue and cultured cells of human and mouse.
Usage	For research use only (RUO)
Storage	Store DB3 and DB5 at –20°C away from light. Store DB4, DB6 and 8 well assay strips at 4°C away from light.Store other components at room temperature. All components are stable for 6 months from date of shipment in proper storage.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Kit Components

DB1 (10X wash buffer) 25 ml DB2 (assay buffer) 3 ml DB3 (Dnmt3B standard, 20 μ g/ml)* 20 μ l DB4 (affinity antibody 100 μ g/ml)* 50 μ l DB5 (detection antibody 200 μ g/ml)* 20 μ l DB6 (developing solution) 12 ml DB7 (stop solution) 6 ml 8 well substrate-coated strip (with frame) 12* For maximum recovery of the products, centrifuge the original vial after thawing prior to opening the cap.

Features & Benefits

Extremely fast procedure, which can be completed within 3 hours. Innovative colorimetric assay, which quantitatively measures Dnmt3B amount with no need for electrophoresis. Strip microplate format makes the assay flexible: manual or high throughput analysis. Simple, reliable, and consistent assay conditions.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA