

Europium Calibration Kit

Cat. No. Kit-0330 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Europium Calibration Kit is used for optimizing and calibrating the europium signal output in all time-resolved fluorescence readers.
Usage	For research use only (RUO)
Storage	Store at 4°C.
Kit Components	Enhancement Solution 50 ml Europium chloride 100 nM

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA