

Fructose-6-Phosphate Assay Kit

Cat. No. Kit-0331 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The PicoProbe F6P Assay Kit can detect F6P in the range of 0.01 to 0.5 nmoles with detection sensitivity ~1 μ M of F6P. (100 assays) Simple procedure; takes ~ 35-40 minutes to complete the assayFast and convenientSensitive, efficient and accurate assay
Size	100 assays
Applications	The PicoProbe F6P Assay Kit can detect F6P in the range of 0.01 to 0.5 nmoles with detection sensitivity ~1 μ M of F6P.
Target Species	Mammalian
Usage	This product is furnished for LABORATORY RESEARCH USE ONLY. Not for diagnostic or therapeutic use.
Kit Components	F6P Assay BufferF6P ProbeF6P Enzyme MixF6P ConverterF6P DeveloperF6PStandard (10 μ mol)
Detection method	Fluorescence Plate ReaderFluorometer
Compatible Sample Types	Cell Culture SupernatantsFood ProductsGrowth MediaPlasmaSerumUrine

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA