

Factor IXa Activity Fluorometric Assay Kit

Cat. No. Kit-0335 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Rapid and convenient The assay can detect Factor IXa activity as low as 10 pg.
Description	The coagulation Factor IX (or Christmas factor) (EC 3.4.21.22) is a vitamin K-dependent serine protease. Factor IX is produced as an inactive precursor and is activated via cleavage by either factor XIa (contact pathway) or factor VIIa (tissue factor pathway). In the presence of calcium ions and negatively charged membrane phospholipids, activated factor IX (FIXa) then binds to the activated Factor VIII (FVIIIa) and proteolytically activates factor X (FX) to factor Xa (FXa). Factor IXa Activity Assay kit is based on the ability of FIXa to generate FXa. The generated FXa proteolytically cleaves a synthetic substrate and releases a fluorophore, AMC, which can be easily quantified by fluorescence microplate reader. The assay is simple, rapid and can detect activity as low as 10 pg of FIXa in a variety of samples.
Applications	Determine activity of purified Factor IXa-Detect the activity of Factor IXa in plasma and serum samples
Usage	For Research Use Only! Not For Use in Humans.
Storage	-20°C
Size	100 assays
Kit Components	• FIXa Assay Buffer • FXa Substrate-AMC

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

- Enzyme Mix I
- Enzyme Mix II
- Phospholipids
- FIXa Enzyme Standard (10 ng)

Detection method	Fluorescence (Ex/Em = 360/450 nm)
-------------------------	------------------------------------

Compatible Sample Types	Purified Factor Ixa, Plasma, Serum
--------------------------------	------------------------------------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA