

GAA Activity Colorimetric Assay Kit

Cat. No. Kit-0353 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	GAA Activity Assay Kit (Colorimetric) is used for measuring α -Glucosidase activity using colorimetric methods.
Description	α -Glucosidase breaks down α -1,4 linked polysaccharides to glucose, which can be utilized as a source of energy. In the biotechnology industry, α -glucosidase is used to produce glucose from intermediate breakdown products of starch hydrolysis generated by enzymes such as amylase. Pompe disease, one of the 12 known glycogen storage diseases, is an autosomal recessive metabolic disorder attributed to α -glucosidase deficiency. In this disease, glycogen accumulates in the lysosomes, resulting in progressive muscle weakness, heart failure and other neurological symptoms. In GAA Activity Assay Kit (Colorimetric), α -Glucosidase hydrolyzes the Substrate Mix to release the p-nitrophenol that can be measured colorimetrically (OD = 410 nm). This is an easy, quick and high-throughput capable kit that can measure 0.1-10 mU of α -glucosidase activity in a variety of samples.
Applications	Measurement of α -Glucosidase activity in biological samples Screening α -Glucosidase inhibitors
Usage	For research use only (RUO)
Storage	Store kit at -20°C, protected from light. Warm α -Glucosidase Assay Buffer to room temperature before use. Briefly centrifuge small vials prior to opening.
Kit Components	α -Glucosidase Assay Buffer 25 ml α -Glucosidase Substrate Mix 300 μ l α -Glucosidase

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Positive Control 1 vialp-Nitrophenol Standard (100 mM) 100 µ	
Detection method	Colorimetric
Compatible Sample Types	Cell culture, Saliva, Serum, Tissue

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA