

Gel Staining/Destaining Kit

Cat. No. Kit-0361 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Gel Staining/Destaining Kit is used for SDS-PAGE gel staining and destaining.
Description	SDS-PAGE is one of the most powerful and commonly used techniques in molecular biology. However, traditional protein gel staining and destaining methods take several hours to finish and leaves a high blue background on the gel reducing the detection sensitivity. Gel Staining/Destaining Kit provides a fast and sensitive method for SDS-PAGE gel staining and destaining. The procedure requires only 5 min for gel staining and 15 min for destaining. A total of 20 min allow visualization of as little as 5 ng of protein on a water-clear background. The kit contains sufficient solutions to stain and destain up to 40 pieces of protein mini gels.
Usage	For research use only (RUO)
Storage	Store Solution A and Solution B at room temperature. Both solutions are stable for at least 1 year.
Kit Components	2X Solution A. Add 250 ml of isopropanol (not provided) before use to generate the 500 ml of 2X Solution A. 500 ml 6X Solution B 500 ml

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA