

HAT Activity Fluorometric Assay Kit

Cat. No. Kit-0408 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Simple, rapid & convenient The assay can detect HAT Activity as low as 6 mU in a variety of samples
Description	Histone Acetyltransferases (HATs) are enzymes that acetylate histone substrates resulting in important regulatory effects on chromatin structure and assembly, and gene transcription. Modifications of these proteins by HATs play an important role in the control of gene expression, and their dysregulation has been linked to cancer, neurodegeneration, and other diseases. HAT Activity Assay Kit utilizes Acetyl CoA and H3 histone peptide as substrates. In this assay, HAT enzyme catalyzes the transfer of acetyl groups from Acetyl-CoA to the histone peptide, thereby generating two products-acetylated peptide and CoA-SH. The CoA-SH reacts with the developer to generate a product that is detected fluorometrically at Ex/Em = 535/587 nm. The assay can detect HAT activity as low as 6 mU in a variety of samples.
Applications	Measurement of HAT activity in Nuclear Extracts Measurement of HAT activity of purified enzyme preparations
Usage	For Research Use Only! Not For Use in Humans
Storage	-80°C
Size	100 assays
Kit Components	• HAT Assay Buffer • Acetyl CoA (Lyophilized)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

- H3 Peptide (Lyophilized)
- Substrate Mix (Lyophilized)
- Developer
- PicoProbe
- CoA Standard (Lyophilized)
- Positive Control (HeLa Nuclear Extract)

Target Species	N/A
-----------------------	-----

Detection method	Fluorescence (Ex/Em = 535/587 nm)
-------------------------	-----------------------------------

Compatible Sample Types	• Nuclear extracts from cells and tissue • Recombinant enzyme
--------------------------------	---

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA