

HDAC3 Inhibitor Screening Kit

Cat. No. Kit-0432 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	HDAC3 Inhibitor Screening Kit is used for screening HDAC3 inhibitor.
Description	Histone deacetylases (HDACs) represent a large family of enzymes identified as key regulators of nucleosomal histone acetylation, a major event that controls eukaryotic gene transcription and are classified into three groups. They are believed to be involved in important biological activities, such as cell differentiation, proliferation, apoptosis, and senescence. In HDAC3 Inhibitor Screening Kit, HDAC-3 and developer will deacetylate and cleave the substrate [Arg-His-Lys-Lys(Ac)-AFC] to release the AFC molecule, which can be detected by fluorometric method (Ex/Em = 380/500 nm). In presence of a HDAC-3 inhibitor, the cleavage will be inhibited. The kit provides a rapid, simple, sensitive, and reliable test, which is also suitable for high throughput screening of HDAC-3 inhibitors. Trichostatin A (TSA) is included as a control inhibitor to compare the efficacy of the test inhibitors.
Usage	For research use only (RUO)
Storage	Store the kit at -20°C, protect from light. Store HDAC-3 enzyme at -80°C. Allow Assay Buffer to warm to room temperature before use. Briefly centrifuge vials before opening. Read the entire protocol before performing the assay.
Kit Components	HDAC3 Assay Buffer. Cap color: WM. 25 ml HDAC3 Substrate. Cap color: red. 200 μl HDAC-3 Enzyme. Cap color: green. 200 μl Developer. Cap color: orange. 1 ml Trichostatin A. Cap color: blue. 50 μl

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Detection method	Fluorometric
-------------------------	--------------