

# Progesterone Receptor Competitor Assay Kit, Green

Cat. No. Kit-0969    Lot. No. (See product label)

## SPECIFICATION

### Product Overview

The Progesterone Receptor Competitor Assay Kit, Green, provides a sensitive and efficient method for high-throughput, fluorescence polarization-based screening of potential progesterone receptor (PR) ligands. The kit uses insect cell-expressed, human PR ligand binding domain tagged with glutathione-S-transferase (GST) (also available separately) and a novel, high-affinity, fluorescent progesterone receptor ligand (PL Green) in a homogenous, mix-and-read assay format.

This kit produces a green color read-out; a "red" kit is also available. Both kits suitably determine the IC<sub>50</sub>s of competitive compounds. Color choice depends on the spectral properties of the compounds to be screened and the available instrumentation.

### Description

The PolarScreen Progesterone Receptor Competitor Assay is:

- Complete & ready-to-use out of the box—all reagents included, just add your test compounds
- High throughput compatible—fluorescence polarization assays in a 384-well plate
- Convenient—components available in bulk sizes for larger screening needs

#### How the Assay Works

An N-terminal fusion of glutathione transferase to the ligand-binding domain of the human progesterone receptor (PR-LBD(GST)) is added to a fluorescently-tagged progesterone ligand (Fluormone PL Green) in the presence of competitor test compounds in microwell plates. The presence of effective competitors prevents the formation of a PL Green/PR-LBD(GST) complex, resulting in a decrease of the polarization value. The shift in polarization value in the presence of test compounds is

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

used to determine relative affinity of test compounds for PR-LBD(GST).  
Contact us for bulk quantities or for assistance with other nuclear receptor assay needs, from assay development to library screening and beyond.

**Applications**

Competitive Binding Assay, Nuclear Receptor Assay

**Storage**

1 tube Progesterone Receptor Ligand Binding Domain Recombinant Protein (store at -68 to -85°C)  
1 tube of Fluormone PL Green (store at -5 to -30°C)  
2 bottles of PR Screening Buffer (store at -5 to -30°C)  
1 tube of 1 M DTT (store at -5 to -30°C)

**Size**

400 x 20 µ assays

 Tel: 1-631-559-9269 1-516-512-3133 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA