

Invertase Activity Colorimetric Assay Kit

Cat. No. Kit-0478 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Invertase Activity Assay Kit (Colorimetric) is used to measure Invertase activity.
Description	Invertase (EC 3.2.1.26), also known as sucrase or β -fructofuranosidase, catalyzes the hydrolysis of sucrose (table sugar) by cleaving its glycosidic bond and forming one molecule each of glucose and fructose. Invertase is widely expressed among plants and microorganisms including yeast and bacteria. It is extensively utilized in industrial processes, including food, ethanol production, cosmetics, drugs, paper and bioelectronics electrodes. The Invertase Activity Assay Kit (Colorimetric) provides a simple, sensitive and high-throughput adaptable assay that detects relevant concentrations of Invertase in biological samples or purified/crude enzyme preparations. The assay is based on the enzymatic detection of glucose, one of the hydrolysis products. Detection limit: 20 μ U.
Applications	Measurement of Invertase activity in various biological samples/preparations.
Target Species	Bacteria, Mammals
Usage	For research use only (RUO)
Storage	Store kit at -20°C, protected from light.
Kit Components	Invertase Hydrolysis Buffer 25 mL Invertase Assay Buffer 25 mL Invertase Stop Solution 1 mL OxiRed Probe (in DMSO) 0.2 mL Invertase Enzyme Mix (Lyophilized) 1 vial Invertase Substrate 1.2 mL Glucose Standard (100 mM) 100 μ L Positive Control

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

	(Lyophilized) 1 vial
Detection method	Colorimetric
Compatible Sample Types	Biological fluids: e.g. SalivaBacteriaCrude/purified enzyme preparations

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA