

Maleimide Green Fluorescence Assay Kit

Cat. No. Kit-0534 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Maleimide Assay Kit (Green Fluorescence) is a one-step fluorometric method for the quantification of maleimide groups in a variety of molecules such as proteins.
Description	Sensitive assays of maleimide and thiol groups are required for the efficient conjugation of proteins that are expensive and available only in small amounts. A variety of crosslinking reagents with a maleimide group are widely used for crosslinking proteins to proteins or proteins to other biomolecules. There are few reagents or assay kits available for quantifying the number of maleimide groups that are introduced into the first protein. All the commercial kits have tedious protocols.
Applications	Maleimide Assay Kit (Green Fluorescence) uses a proprietary dye that has enhanced fluorescence upon reacting with a maleimide.
Usage	For research use only (RUO)
Storage	Store the kit at -20°C and avoid from light.
Kit Components	Component A: Maleimide Green 1 vial Component B: Reaction Buffer 500 µL Component C: Assay Buffer 25 mL Component D: N-ethylmaleimide Standard (10 mM) 50 µL Component E: DMSO 200 µL
Detection method	Fluorometric
Features & Benefits	Broad Application: Can be used for quantifying maleimide group in a variety of molecules such as proteins. Sensitivity: Detects as low as 10 picomoles of

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

maleimide. Continuous: Easily adapted to automation without a separation step. Convenient: Formulated to have minimal hands-on time. No wash is required. Non-Radioactive: No special requirements for waste treatment.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA