

Maleimide Rapid Assay Kit

Cat. No. Kit-0535 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Maleimide Assay Kit (Rapid) is a quantitative assay for the quantification of maleimide groups.
Description	A variety of crosslinking reagents with a maleimide group are widely used for crosslinking proteins to proteins, or proteins to other biomolecules. A key challenge in maleimide crosslinking technology has been the development of reliable and reproducible methods for quantitation the amount of maleimide linked to a single protein or antibody. Maleimides can be directly assayed spectrophotometrically at 302 nm. However, the small extinction coefficient of 620 M-1cm-1 renders this assay insensitive, and the assay is further complicated by the protein absorbance at the same wavelength. There are a few methods available for quantitating the number of maleimide groups by indirect ways. And in many cases, the maleimide ratio is only estimated by experience.
Applications	Maleimide Assay Kit (Rapid) provides an accurate method to quantify maleimide group by our maleimide sensor Maleimide Blue with the maximum absorbance at ~780nm.
Usage	For research use only (RUO)
Storage	Upon receipt, store Maleimide Blue (Component A) the kit at 2-6°C, kept from light and moisture. When stored properly, the kit components should be stable for six months. Note 1: Do not freeze Spin Column (Component C). Note 2: Warm all the components before run the required assays. 50 to 100 µg maleimide-linked antibody

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

or protein sample is needed for determining the amount of maleimide.

Kit Components

Component A: Maleimide Blue (One vial is for 50 ~100 µg maleimide-linked antibody or protein) 2 vials
Component B: Assay Buffer 15 mL
Component C: Spin Column 2 columns
Component D: Washing Tube 2 mL
Component E: Collecting Tube 1.5 mL

Detection method

Colorimetric

Compatible Sample Types

Polymers, Proteins

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA