

MMP-3 Inhibitor Screening Kit

Cat. No. Kit-0590 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The kit can be used to screen inhibitors of MMP-3. (100 assays) Simple procedure; takes ~ 35-45 minutesFast and convenientThe assay is simple, reliable and ensures accuracy and consistency when working with multiple samples.
Size	100 assays
Usage	This product is furnished for LABORATORY RESEARCH USE ONLY. Not for diagnostic or therapeutic use.
Kit Components	MMP-3 Assay BufferMMP-3 SubstrateMMP-3 Enzyme (lyophilized)Inhibitor Control (0.1 mM)
Detection method	Fluorescence Plate ReaderFluorometer

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA