

Nuclear Extraction Kit

Cat. No. Kit-0635 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Nuclear Extraction Kit is a method for preparing nuclear extraction.
Description	Preparation of nuclear extracts is the first step in examining transcription factor activity. This Nuclear Extraction Kit is formulated for the quick and simple isolation of nuclear and cytoplasmic fractions from cultured cells and tissue homogenates. The proteins isolated using this kit can also be used in electrophoretic mobility shift assays (EMSA) and western blotting applications.
Usage	For research use only (RUO)
Storage	This kit will perform as specified if stored as directed and used before the expiration date indicated on the outside of the box. Item Storage Nuclear Extraction Hypotonic Buffer (10X) 4°C Nuclear Extraction Dithiothreitol (1 M) -20°C Nuclear Extraction Protease Inhibitor Cocktail (100X) -20°C Nuclear Extraction PBS (10X) 4°C Nuclear Extraction Phosphatase Inhibitors (50X) -80°C Nuclear Extraction Buffer (2X) 4°C Nonidet P-40 Assay Reagent (10%) RT
Kit Components	Nuclear Extraction Hypotonic Buffer (10X) 6 mL Nuclear Extraction Dithiothreitol (1 M) 60 µL Nuclear Extraction Protease Inhibitor Cocktail (100X) 650 µL Nuclear Extraction PBS (10X) 100 mL Nuclear Extraction Phosphatase Inhibitors (50X) 5 mL x 2 Nuclear Extraction Buffer (2X) 3 mL Nonidet P-40 Assay Reagent (10%) 6 mL

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA