

Deproteinizing Sample Preparation Kit

Cat. No. Kit-0726 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Deproteinizing Sample Preparation Kit is used for removing proteins from samples.
Description	The presence of protein and various enzyme activities frequently interferes with the analysis of small molecules in biological samples. Many bioassays require removal of proteins from samples prior to analysis. Among the deproteinization protocols developed over the last half century, perchloric acid (PCA) precipitation has been extensively used in many different sample preparation procedures, since not only does it remove most of the protein present but it also functions to stabilize many of the small molecule analytes. PCA deproteinization methods have been successfully used in the preparation samples prior to quantitation of an array of small molecules, including glycogen, ATP, cAMP, glutathione, antioxidants, etc. The newly designed Deproteinizing Sample Preparation Kit utilizes a PCA precipitation method, providing a unique tool to prepare samples for the analysis of a variety of small molecules. Using the kit, proteins are precipitated, excess PCA is removed, samples are neutralized. Samples prepared using this kit can be directly used in a wide variety of bioassays. The method is easy, convenient, and can be used for the preparation of a large number of samples.
Target Species	Mammals
Usage	For research use only (RUO)
Storage	Store kit at 4°C. Keep kit components ice while in use. There may be some precipitation in the neutralization buffer. Shake bottle gently a few times to resuspend

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

before using.

Kit Components	PCA. Cap code: WM. 20 mlNeutralization Solution. Cap code: NM. 4 ml
Compatible Sample Types	Cell lysate, Serum, Tissue Homogenate, Urine

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA