

TACE Inhibitor Screening Kit

Cat. No. Kit-0810 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	In TACE inhibitor screening Kit, TACE (tumor necrosis factor- α -converting enzyme) hydrolyzes the specific FRET substrate to release the quenched fluorescent group, which can be detected fluorometrically at Ex/Em318/449 nm. In the presence of the potent TACE inhibitor, the hydrolyzation of substrate will be impeded. The kit provides a rapid, simple, sensitive and reliable test suitable for high-throughput screening of TACE inhibitors and can be modified to check the relative TACE activity. Inhibitor Control GM6001 is included to compare the efficacy of test inhibitors.
Size	100 tests
Description	The TACE (tumor necrosis factor- α -converting enzyme), also called ADAM metalloproteinase domain 17 (ADAM17), is a 70-kDa enzyme that belongs to the ADAM protein family of disintegrins and metalloproteases. TACE is believed to be involved in the processing of tumor necrosis factor alpha (TNF- α) at the surface of the cell, and from within the intracellular membranes of the trans-Golgi network. This process, which is also known as shedding, involves the cleavage and release of a soluble ectodomain from membrane-bound pro-proteins (such as pro-TNF- α), and is of known physiological importance.
Applications	Functional Studies more details
Storage	Store at -20°C.
Kit Components	Components: Identifier; TACE Assay Buffer: WM; TACE Enzyme: Green; TACE

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Inhibitor Control: Purple; TACE Substrate: Red

Assay time	0h 20m
-------------------	--------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA