

TBARS Assay kit II

Cat. No. Kit-0812 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	TBARS Assay kit is used for the quantitative determination of thiobarbituric acid reactive substances through the detection of an adduct of malondialdehyde (MDA) and thiobarbituric acid (TBA).
Description	The sensitivity of measuring Thiobarbituric Acid Reactive Substances (TBARS) has made this assay the method of choice for screening and monitoring lipid peroxidation, a major indicator of oxidative stress. Oxidative stress in the cellular environment results in the formation of highly reactive and unstable lipid hydroperoxides. Decomposition of the unstable peroxides derived from polyunsaturated fatty acids results in the formation of malondialdehyde (MDA), which can be quantified colorimetrically following its controlled reaction with thiobarbituric acid. This assay remains the most widely employed assay used to determine lipid peroxidation.
Applications	The TBARS Assay kit measures malondialdehyde (MDA), a reactive compound formed from lipid peroxides that are generated under conditions of oxidative stress. Oxidative modification of lipids occurs with aging and various diseases, and increased oxidative stress is associated with diabetes and its complications. MDA forms an adduct with thiobarbituric acid (TBA). Results are calculated from a standard curve constructed with authentic MDA. Exocells TBARS assay can be used with a spectrum of biological samples including body fluids, tissue and cell specimens.
Usage	For research use only (RUO)
Storage	Store all kit reagents at 2-8°C. The components should be used before the

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

	expiration date indicated on the outside of the box. TBA stock should be refrigerated.
Kit Components	Microplate for fluorometer. 96 (8x12) wells Thiobarbituric Acid stock solution 1 vial TBARS Assay Diluent 1 vial MDA Standard 10 μ M Malondialdehyde 1 vial Bis(dimethylacetal) 1 vial Acetic Acid 1 vial
Detection method	Fluorimetric/Colorimetric
Compatible Sample Types	Body fluids, Cell specimens and Tissue

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA