

Beta-Lactamase Activity Colorimetric Assay Kit

Cat. No. Kit-0861 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Beta-Lactamase Activity Assay Kit (Colorimetric) is used to detect and quantify the beta-lactamase activity at 490 nm.
Description	Beta-Lactamases (βLS) are a large family of hydrolases comprising more than 850 identified members expressed in Gram-positive and Gram-negative bacteria. βLS can be classified according to their substrate or inhibitor specificity. These enzymes are capable of hydrolyzing four atom rings known as β-lactams. Antibiotics containing β-lactam rings (i.e. penicillin, cephalosporin, monobactam, carbapenem) are highly susceptible to be hydrolyzed via enzymatic activity, which deactivates their antibiotic potency. βLS have become a significant clinical threat due to the alarming number of cases of bacterial strains showing β-lactam antibiotic resistance. Beta-Lactamase Activity Assay Kit offers a simple and sensitive assay that can detect and quantify the enzymatic activity of these hydrolases. The assay is based on the hydrolysis of Nitrocefin, a chromogenic cephalosporin, that results in the generation of a colored product (OD 490 nm), which is directly proportional to the amount of βL activity. The assay can detect enzymatic activity as low as 0.06 mU in a variety of biological samples.
Applications	Measurement of β -Lactamase activity in various biological samples Analysis of β -Lactamase activity in pathological conditions
Usage	For research use only (RUO)
Storage	Store the kit at -20°C, protected from light.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Kit Components	β L Assay Buffer 27 mL Nitrocefin (in DMSO) 220 μ Positive Control (Lyophilized) 1 vial β L Hydrolysis Buffer 100 μ
Detection method	Colorimetric
Compatible Sample Types	Bacterial cultures, Fermentation media, etc., Food e.g. milk, Saliva, Serum, Urine from mammals infected with β -secretase secreting bacteria

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA