

DNA Laddering Kit

Cat. No. Kit-1945 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Size	24 reactions
Description	Apoptosis is associated with the fragmentation of chromosomal DNA into multiples of the 180 bp nucleosomal unit, known as DNA laddering. This kit is designed to maximize extraction of small fragments of DNA with minimal contamination from intact chromatin. The isolated DNA is separated by electrophoresis and can be visualized using ethidium bromide.
Storage	-20°C
Kit Components	Lysis Buffer: 4 vials/1.2 ml; -20°C 10% SDS: 1 vial/0.48 ml; -20°C Enzyme A: 1 vial/0.48 ml; -20°C Enzyme B: 1 vial/0.48 ml; -20°C Precipitating Reagent: 3 vials/1.04 ml; -20°C Gel-Loading Buffer (6X): 1 vials/0.48 ml; -20°C

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA