

Nitric Oxide Fluorimetric Detection Kit (Red)

Cat. No. Kit-2046 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Nitrixyte probes are developed and used in our kit as an excellent replacement for DAF-2 for the detection and imaging of free NO in cells. Compared to the commonly used DAF-2 probe, Nitrixyte probes have better photo stability and enhanced cell permeability. This particular kit uses Nitrixyte Red that can react with NO to generate a bright red fluorescent product that has spectral properties similar to Texas Red. Nitrixyte Red can be readily loaded into live cells, and its fluorescence signal can be conveniently monitored using the filter set of Red. This kit is optimized for flow cytometry applications.
Size	100 assays
Description	Nitric oxide (NO) is an important biological regulator involved in numbers of physiological and pathological processes. Altered NO production is implicated in various immunological, cardiovascular, neurodegenerative and inflammatory diseases. As a free radical, NO is rapidly oxidized and there is relatively low concentrations of NO existing in vivo. It has been challenging to detect and understand the role of NO in biological systems. Fluorimetric Intracellular Nitric Oxide Assay Kit provides a sensitive tool to monitor intracellular NO level in live cells.
Applications	The kit is optimized for flow cytometry applications.
Storage	Keep Component A&B at -20°C and Component C at 4°C.
Kit Components	Component A: Nitrixyte Red (500x) 100 μ l Component B: NONOate Positive Control 1

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

vial (lyophilized powder) Component C: Assay Buffer 1 bottle (10 ml)

Detection method	Flow cytometer
-------------------------	----------------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA