

Bile Acid Assay Kit

Cat. No. Kit-2070 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Bile Acid Assay Kit provides a convenient fluorimetric means to measure total bile acids in biological samples. In the assay, 3 α -hydroxysteroid dehydrogenase reacts with all twelve bile acids, converting NAD to NADH, which reduces a probe to a highly fluorescent product. The resulting fluorescence intensity ($\lambda_{exc/em}$ = 530/585 nm) is linear to the bile acid concentration in the sample.
Size	100 tests
Description	Twelve different types of bile acids are typically found in mammals, among them two primary types are cholic acid and chenodeoxycholic acid. These can be dehydroxylated into secondary bile acids. Finally, these four can be conjugated to either taurine or glycine creating 8 different conjugated bile acids. Bile acid levels in feces, blood, urine, and bile can be used as markers for various diseases such as hyperlipidemia, cholestasis, gall stones, colon cancer, etc. Bile acids also exist in a sulfate salt form known as bile acid-sulfates. Sulfation of bile acids increases their solubility and decreases intestinal absorption, thereby enhancing fecal and urinary excretion. This assay does not measure bile acid-sulfates, and measures only the twelve non-sulfated bile acids.
Applications	For quantitative determination of total bile acids and evaluation of drug effects on bile acid metabolism.
Storage	The kit is shipped on ice. Store all reagents at -20°C. Shelf life of 6 months after receipt.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Kit Components	Assay Buffer: 10 mL Enzyme A: 120 μ NAD Solution: 1 mL Enzyme B: 120 μ Probe: 750 μ Standard: 120 μ
Detection method	FL530/585nm
Compatible Sample Types	Serum, plasma, urine and other biological samples.
Features & Benefits	• Safe. Non-radioactive assay. • Sensitive and accurate. Linear detection range of 1 - 150 μM bile acids. • Convenient and high-throughput. Homogeneous "mix-incubate-measure" type assay. No wash and reagent transfer steps are involved. Can be readily automated on HTS liquid handling systems for processing thousands of samples per day.
Assay time	20 min
Sensitivity	Detection Limit: 1 μ M

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA