

α-Mannosidase Assay Kit

Cat. No. Kit-2079 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The nonradioactive, colorimetric AMA assay is based on the cleavage of 4-nitrophenol from the synthetic substrate. Nitrophenol becomes intensely colored after addition of the stop reagent. The increase in absorbance at 405 nm after addition of the stop reagent is directly proportional to the enzyme activity.
Size	100 tests
Description	α-MANNOSIDASE (AMA) is an enzyme which catalyzes the cleavage of the alpha form of mannose. α-Mannosidase assists in the breakdown of complex sugars from glycoproteins in the lysosome. Defective AMA or deficient AMA activity causes α-mannosidosis and leads to deterioration of the central nervous system in children.
Applications	For quantitative determination of α-mannosidase enzyme activity using a p-nitrophenol based substrate.
Storage	The kit is shipped at room temperature. Store all components at 4°C upon receiving. Shelf life: 6 months after receipt.
Kit Components	Substrate: 10 mL Stop Reagent: 12 mL Standard: 1 mL 12.5 mM Nitrophenol
Detection method	OD405nm
Compatible Sample Types	Biological samples e.g. plasma, serum, tissue and culture media.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Features & Benefits	• Fast and sensitive. Linear detection range (10 μL sample): 1 to 250 U/L for a 10 minute reaction. • Convenient and high-throughput. Homogeneous "mix-incubate-measure" type assay. Can be readily automated on HTS liquid handling systems for processing thousands of samples per day.
--------------------------------	--

Assay time	10 min
-------------------	--------

Sensitivity	Detection Limit: 1 U/L
--------------------	------------------------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA