

Xanthine Assay Kit

Cat. No. Kit-2081 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Simple, direct and high-throughput assays for measuring xanthine find wide applications in research and drug discovery. The xanthine assay kit uses a single Working Reagent that combines the xanthine oxidase reaction and color reaction in one step. The change in color intensity of the reaction product at 570 nm or fluorescence intensity at $\lambda_{ex}/\lambda_{em}$ = 530/585 nm is directly proportional to xanthine concentration in the sample.
Size	100 tests
Description	Xanthine is a purine base that can be found in most animal tissues and fluids. It is a product in the purine degradation pathway, produced by guanine deaminase from guanine, and by xanthine oxidoreductase from hypoxanthine. Xanthine is degraded to uric acid by xanthine oxidase. Clinically, xanthine and its derivatives act on sleep-inducing adenosine receptors as antagonists.
Applications	For quantitative determination of xanthine and drug effects on its metabolism.
Storage	The kit is shipped on ice. Store all reagents at -20°C. Shelf life: 6 months after receipt.
Kit Components	Assay Buffer: 10 mL HRP Enzyme: 120 μ Standard: 1 mL 2 mM Xanthine Dye Reagent: 120 μ XO Enzyme: 100 μ
Detection method	OD570nm, FL530/585nm
Compatible Sample	Cell lysate, serum, and other biological samples.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Types**Features & Benefits**

• Sensitive and accurate. Use as little as 10 μ samples. Linear detection range in 96-well plate for 30 minute incubation: 0.01 to 2.5 mM xanthine for colorimetric assays and 3 to 250 μ M for fluorimetric assays. • Simple and convenient. The procedure involves addition of a single working reagent and incubation for 30 min at room temperature. • Fast and high-throughput. Assays using 96-well plates and liquid handling system could allow simultaneous processing tens of thousands of samples per day.

Assay time

30 min

SensitivityDetection Limit: FL 3 μ M; OD 10 μ M Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA