

Acidic Mammalian Chitinase Inhibitor Screening Kit

Cat. No. Kit-2104 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Acidic Mammalian Chitinase Inhibitor Screening Kit uses Bisdionin F— a highly selective AMCase inhibitor as positive control and it provides a simple, rapid, sensitive and reliable test suitable for screening/characterizing AMCase inhibitors.
Size	100 assays
Description	Acidic Mammalian Chitinase (AMCase) belongs to the glycosyl hydrolase-18 family (EC 3.2.1.14), which is capable of hydrolyzing the β - (1, 4)-linkage between the adjacent N-acetyl glucosamine residues of chitin. AMCase is relatively abundant in the gastrointestinal tract and lung. It is highly expressed in the lungs of asthmatic patients and asthma animal models. In addition, the change in AMCase expression can be used to monitor Th-2/IL-13 mediated immune response and inflammation. Inhibition of AMCase activity has been considered as a therapeutic strategy to alleviate symptoms of allergic airway inflammation. Thus, there is a need to identify AMCase inhibitors as potential candidates to treat and/or prevent asthma.
Applications	Screening/characterizing/studying Mammalian Chitinase inhibitors/activators
Storage	Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay. AMCase Assay Buffer: Warm to 37 °C before use. Store at either 4°C or -20°C. AMCase Substrate: Store at -20°C. Bring to room temperature before use. Acidic Mammalian Chitinase: Reconstitute Acidic Mammalian Chitinase in 550 μ l AMCase Assay Buffer and mix thoroughly. Aliquot and store at -20°C. Avoid repeated freeze/thaw. Keep on ice while in use. Use within

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

two months. Bisdionin F: Store at -20°C. Bring to room temperature before use.

Kit Components

AMCase Assay Buffer: 25 ml
AMCase Substrate (in DMSO): 25 µl
Acidic Mammalian Chitinase (lyophilized): 1 vial
Bisdionin F (5 mM): 20 µl

Detection method

Fluorescence (RFU= 320/445 nm)

Features & Benefits

- Simple and reliable test to screen AMCase inhibitors/activators
- High-throughput compatible
- Includes Inhibitor Control (Bisdionin F)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA