

HRV 3C Protease Inhibitor Screening Kit

Cat. No. Kit-2118 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	HRV 3C Protease Inhibitor Screening Kit utilizes the ability of a 3C Protease (derived from a HRV rhinovirus-14, EC: 3.4.22.28) to cleave a chromogenic peptide substrate to release a chromophore (pNA) which can be easily quantified using a microplate reader. In the presence of a HRV 3C Protease-specific inhibitor, the cleavage of the substrate is reduced/abolished resulting in decrease or total loss of the absorbance. This simple and high-throughput adaptable assay kit can be used to screen/study/characterize potential inhibitors of HRV 3C Protease.
-------------------------	---

Size	100 assays
-------------	------------

Description	Human rhinovirus (HRV) infections are the most frequent causative agents of common cold and various other upper respiratory tract infections. Rhinoviruses are members of the picornavirus family, which have a positive-sense, single-stranded RNA genome that is translated into a single polyprotein precursor. In the case of HRVs, the viral polyprotein is mainly processed by the proteases (2A and 3C) to generate functional proteins and enzymes.
--------------------	---

Applications	Screening/studying/characterizing inhibitors of HRV 3C Protease.
---------------------	--

Storage	Store kit at -20°C, protected from light. Briefly centrifuge small vials at low speed prior to opening. Read the entire protocol before performing the experiment. HRV 3C Protease Assay Buffer: Store at 4°C. Bring to room temperature before use. HRV 3C Protease: Aliquot and store at -80°C for long term. Avoid repeated freeze/thaw. Each vial contains enough enzyme solution for 50 assays.
----------------	--

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Kit Components	HRV 3C Protease Assay Buffer: 25 ml HRV 3C Protease: 2 x 100 µl Substrate: 500 µl HRV 3C Protease Inhibitor (2 mM): 20 µl
Detection method	Absorbance (OD= 405 nm)
Compatible Sample Types	Small molecules
Features & Benefits	• Simple, rapid assay. • HTP adaptable.

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA