

Sulfatase Colorimetric Activity Assay Kit

Cat. No. Kit-2146 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Sulfatase Activity Assay Kit provides a quick and easy way to determine sulfatase activity. The kit measures the hydrolysis of a sulfate ester to 4-nitrocatechol, which can be detected at 515 nm. The kit is suitable for measuring activity of purified enzyme as well as sulfatase from biological samples. The limit of detection is below 1 mU.
Size	100 assays
Description	Sulfatases (EC 3.1.6) are enzymes in the esterase class that catalyze the hydrolysis of sulfate esters from a wide range of biological molecules, including steroids, carbohydrates, and proteins. They can be found in intracellular and extracellular spaces and are distributed in a wide range of cells and tissues. Intracellular sulfatases are commonly found localized within the lysosome. Genetic defects in sulfatase can result in certain lysosomal storage disorders and abnormal expression can contribute to certain hormone-dependent cancers, such as breast and prostate cancer.
Applications	Measurement of sulfatase activity
Storage	Store kit at -20°C , protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay. Sulfatase Assay Buffer, Sulfatase Substrate, Stop/Developing Solution, and 4-Nitrocatechol Standard: Bring to room temperature before use. Sulfatase: Reconstitute with 20 $\mu\text{dH}_2\text{O}$. Aliquot and store at -20°C . Avoid repeated freeze/thaw. Stable for two months. Keep on ice during use.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Kit Components	Sulfatase Assay Buffer: 5 mlSulfatase Substrate: 4 mlStop/Developing Solution: 10 mlSulfatase: 1 Vial4-Nitrocatechol Standard (0.5 mM): 1.5 ml
Detection method	Absorbance (515 nm)
Compatible Sample Types	• Animal tissues such as liver, endometrium, and breast cancer• Adherent or suspension cells• Purified enzyme
Features & Benefits	• Simple, Sensitive & Rapid protocol• Highly-Sensitive: Detects as low as 1 mU sulfatase in various sample types

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA