

Cross-linked Allophycocyanin protein

Cat. No. XL-APC-136 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Native Cross-linked Allophycocyanin protein. Absorption Maximum: 650 nm± 2nm. Emission Maximum: 660 nm± 2nm. Extinction Coefficient: 7.0 x 10 ⁵ cm ⁻¹ M ⁻¹ .
Description	Cross-linked Allophycocyanin (XL-APC) is the brightest and most stable form of APC available. The αβ subunits of the APC have been crosslinked in order to provide increased structural integrity in the presence of chaotropic salts.
Form	Supplied in 100 mM sodium phosphate buffer pH7.0 with 60% ammonium sulfate, 1mM EDTA, and 1mM sodium azide.
Molecular Mass	105 kDa
Purity	A650/A280 ≥ 5.0 A650/A620 ≥ 1.5 Cross-link ratio ≥ 1.2
Notes	XL-APC is stored in ammonium sulfate solution. Therefore, XL-APC needs to be re-suspended until no visible precipitates. Firstly, it must be extensively dialyzed prior to use or column to remove any remaining ammonium sulfate. Use PBS or your conjugation buffer (1 ml of XL-APC to 1 liter of buffer) to dialyze at least two changes of buffer.
Stability	2 years under proper storage conditions.
Storage	Store at 4 centigrade, away from light. Do not freeze.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

For research or further manufacturing use only.

Concentration	15 mg/ml
----------------------	----------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA