



PRMT6 (HRMT1L6) Chemiluminescent Assay Kit

Product Information

Cat.No.

Kit-1860

Product Overview

PRMT6 is a histone methyltransferase that exhibits methylation activity toward H3-R2 and H4-R3.

Size

96 reactions

Description

The PRMT6 Direct Activity Assay kit is designed to measure PRMT6 activity for screening and profiling applications. The key to the PRMT6 Direct Activity Assay Kit is a highly specific antibody that recognizes methylated R3 residue of Histone H4. With this kit, only three simple steps on a microtiter plate are required for methyltransferase detection. First, S-adenosylmethionine is incubated with a sample containing assay buffer and methyltransferase enzyme for one hour. Next, primary antibody is added. Finally, the plate is treated with an HRP-labeled secondary antibody followed by addition of the HRP substrate to produce chemiluminescence that can then be measured using a chemiluminescence reader.

Applications

Great for studying enzyme kinetics and screening small molecular inhibitors for drug discovery and HTS applications.

Storage

12 months from date of receipt, when stored as directed. Kit components require different storage conditions. Be sure to store each component at the proper temperature upon arrival.

Kit Components

PRMT6 human recombinant enzyme: 20 µg; -80°C 20 µM S-adenosylmethionine*: 250 µl; -80°C Primary antibody 4: 100 µl; -80°C Secondary HRP-labeled antibody 2: 10 µl; -80°C 4x HMT assay buffer 6: 3 ml; -20°C Blocking buffer: 50 ml; +4°C HRP chemiluminescent substrate A (translucent bottle): 6



CREATIVE **BIOMART**[®]
Assay Kit

PRMT6 (HRMT1L6) Chemiluminescent Assay Kit

ml; +4°C HRP chemiluminescent substrate B (brown bottle): 6 ml; +4°C 96-well plate precoated with histone substrate: 1 plate; +4°C* Decreasing S-adenosylmethionine concentration will make the assay more sensitive to the inhibitors.

Tel: 1-631-559-9269 1-516-512-3133

Email: info@creative-biomart.com

Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA