

## Human Protein C Reference standard

Cat. No. C-19H Lot. No. (See product label)

### SPECIFICATION

<b>Product Overview</b>	This standard is intended for use in the estimation of functional chromogenic activity and antigenic content of protein C concentrates.
<b>Species</b>	Human
<b>Usage</b>	No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution. Allow ampoules to warm to room temperature. Open ampoule, taking care to ensure that all material is in the lower part, and reconstitute with 1.0 ml distilled water. Stand for 10 minutes at room temperature to allow complete dissolution of the material before use. The reconstituted standard should be transferred to a plastic tube and stored on melting ice. The reconstituted Standard should be used as soon as possible, preferably within 4 hours of reconstitution.
<b>Notes</b>	This preparation is not for administration to humans. The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratory's safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.
<b>Stability</b>	Reference materials are held at Creative Biomart within assured, temperature-controlled storage facilities and should be stored on receipt as indicated on the label.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA



Once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use.

**Storage**

Unopened ampoules should be stored in the dark at or below  $-20^{\circ}\text{C}$ . Please note: because of the inherent stability of lyophilized material, Creative Biomart may ship these materials at ambient temperature.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA