

HBsAg Agarose

Cat. No. HBsAg-001A Lot. No. (See product label)

SPECIFICATION

| | |
|------------------------------|--|
| Product Overview | HBsAg Agarose is widely used for the purification of Surface antigen of hepatitis B vaccine. |
| Cat.no | HBsAg-001A |
| Size | 25 mL, 100 mL, 1 L |
| Matrix | Highly cross-linked agarose, 6% |
| Average particle size | 90 µm |
| Ligand | HBsAg |
| Ligand concentration | 10-40 µmol / ml wet gel |
| Recommended flow rate | <150> |
| Viscosity | store at 4-30°C |
| Chemical stability | The quality of resin doesn't change at all when dip in 1 M NaOH, 0.01 M HCl, 6 M guanidine hydrochloride and 8 M urea for 7 days under room temperature. |
| pH working range | 3-13 |

 Tel: 1-631-559-9269 1-516-512-3133


 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA



| | |
|---------------------|------|
| pH CIP range | 2-14 |
|---------------------|------|

| | |
|----------------------------------|---------|
| Temperature Stability | 4-30 °C |
|----------------------------------|---------|

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA