



guanidine carbonate and 10 mM ATP when resuspended. 10X strength. ATP stock: 1x 1 ml, Lyophilized. 100 mM ATP when resuspended. Tris-HCl pH 7.5: 1 x 5 ml, Lyophilized. 100 mM when reconstituted. 96 well plate: 1 plate, Black polystyrene assay plate.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA