

## Ultra Sensitive Nitric Oxide Synthase Assay Kit

Cat. No. Kit-0632 Lot. No. (See product label)

### SPECIFICATION

**Size**

1 kit

**Storage**

1. Store this kit and its components according to the label. The Nitrate Reductase, NADPH Part A and NADPH Part B should be stored at -20°C. The reconstituted Nitrate Reductase may be stored at -20°C for up to six months. The diluted standard curve may be stored at 4°C for up to 24 hours. All other kit components should be stored at 4°C.

**Kit Components**

Component-Contents-Quantity-Storage  
Reaction Buffer: 20 mM HEPES; 0.5 mM EDT A, 50 mL, 4°C; NADPH Part A: NADP<sup>+</sup>, Glucose 6-Phosphate, L-Arginine, 1.0 mL, -20°C; NADPH Part B: Glucose 6-Phosphate dehydrogenase, 1.0 mL, -20°C; Nitrate Reductase: Nitrate Reductase, 1 U, -20°C; Nitrate Reductase Buffer: Proprietary buffer, 1.2 mL, 4°C; Color Reagent 1: Sulfanilamide in 3N HCl, 7.0 mL, 4°C; Color Reagent 2: N-(1-Naphthyl) ethylenediamine dihydrochloride, 7.0 mL, 4°C; Nitrite Standard: 500 pmol/ $\mu$  NaNO<sub>2</sub>, 1.5 mL, 4°C; Microtiter Plates: 96-well low binding, flat-bottom plate, 1, RT

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA