

## PARP1 CoIP Agarose beads

**Cat. No.** PARP1-01A    **Lot. No.** (See product label)

### SPECIFICATION


<b>Product Overview</b>	Immunoprecipitation of PARP1 with PARP1-Trap Agarose. <ul style="list-style-type: none"><li>• Effectively immunoprecipitates PARP1 (human)</li><li>• Functional after IP: on-bead parylation assay</li><li>• No heavy &amp; light chains in your downstream application</li></ul>
<b>Specificity/Target</b>	PARP1, Epitope: Within DNA binding domain of PARP1 does not bind to PARP2, 3, and 9
<b>Tested Reactivity</b>	Human
<b>Conjugate</b>	Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)
<b>Elution buffer</b>	SDS sample buffer 0.2 M glycine pH 2.5
<b>Type</b>	Nanobody
<b>Storage Buffer</b>	20% ethanol
<b>Applications</b>	IP, CoIP

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Western Blot



 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA