

# Recombinant Full Length Panax Ginseng Atp Synthase Subunit A, Chloroplastic(Atpi) Protein, His-Tagged

**Cat. No.** RFL12299PF    **Lot. No.** (See product label)

## SPECIFICATION

<b>Product Overview</b>	Recombinant Full Length Panax ginseng ATP synthase subunit a, chloroplastic(atpl) Protein (Q68S18) (1-247aa), fused to N-terminal His tag, was expressed in E. coli.
<b>Species</b>	Panax ginseng (Korean ginseng)
<b>Source</b>	E.coli
<b>ProteinLength</b>	Full Length (1-247)
<b>Form</b>	Lyophilized powder
<b>AA Sequence</b>	MNVLSCSINTLKGLYDISGVEVGQHFYWKIGGFQVHGQVLITSWVVIAILLGSATI AVR N PQT IPTGGQNFFEYVLEFIRDVSKTQIGEEYGPWVFIGTMFLFIVSNWSGALLP WKII QLPHGELAAPTNDINTTVALALLTSVAYFYAGLTKKGLGYFGKYIQPTPILLPINIL EDF TKPLSLSFRLFGDILADELVVVVLVSLVPSVVPVPIVMFLGLFTSGIQALIFATLAAA YIG ESMEGHH
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Storage</b>	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

**Storage Buffer** Tris/PBS-based buffer, 6% Trehalose, pH 8.0

**Reconstitution** We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## GENE INFORMATION

**Gene Name** atpI

**Synonyms** atpI; PSC0159; ATP synthase subunit a, chloroplastic; ATP synthase F0 sector subunit a; F-ATPase subunit IV

**UniProt ID** [Q68S18](#)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA