

# Recombinant Full Length Danio Rerio Fat Storage-Inducing Transmembrane Protein 2(Fitm2) Protein, His-Tagged

Cat. No. RFL15122DF Lot. No. (See product label)

## SPECIFICATION

<b>Product Overview</b>	Recombinant Full Length Danio rerio Fat storage-inducing transmembrane protein 2(fitm2) Protein (Q52KL1) (1-252aa), fused to N-terminal His tag, was expressed in E. coli.
<b>Species</b>	zebrafish
<b>Source</b>	E.coli
<b>ProteinLength</b>	Full Length (1-252)
<b>Form</b>	Lyophilized powder
<b>AA Sequence</b>	MAAAVAGSLVDKLVCLWRQPYTRIIYLPHLFFCISLVGSLVKNELVPESYFSSSRNVL NL YFVKVSWGWTIVLLLPIAYSNFYIKSHMFALRRLTSLLVATLVWYICTETFFYIEDI TG SCYESNTMVVIRGEFDTKAACRKAGFFWDGFDISGHSFILSYSSLVIMEEMVPML HIQPA YRNPLDCLYLALNVIVAIWIWMFGCTSVYFHDIIDKILGTSCGILGWYMTYKV WYVKLF SPGLPPQPKQHT
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Storage</b>	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

repeated freeze-thaw cycles.

**Storage Buffer** Tris/PBS-based buffer, 6% Trehalose, pH 8.0

**Reconstitution** We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## GENE INFORMATION

**Gene Name** fitm2

**Synonyms** fitm2; fit2; zgc:110840; Acyl-coenzyme A diphosphatase FITM2; Fat storage-inducing transmembrane protein 2; Fat-inducing protein 2

**UniProt ID** [Q52KL1](#)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA