

Recombinant Full Length Uukuniemi Virus Envelope Glycoprotein(Gp) Protein, His-Tagged

Cat. No. RFL23709UF Lot. No. (See product label)

SPECIFICATION

Product Overview	Recombinant Full Length Uukuniemi virus Envelope glycoprotein(GP) Protein (P09613) (514-1008aa), fused to N-terminal His tag, was expressed in E. coli.
Species	Uukuniemi virus
Source	E.coli
ProteinLength	Full Length of Mature Protein (514-1008)
Form	Lyophilized powder
AA Sequence	<p>CSDSISVTASSQRCSTSSDGVNSCFVSTSSLLQVSPKGQESCLILKGPTGTAVDSIRI KT TDIKLECVRRDLYWVPRVTHRCIGTRRCHLMGACKGEACSEFKINDYSPEWGHE EELMAQ LGWSYCVEQC GGALCQCFNMRPSCFYLRKTF SHLSQDAFNIYECSEWS YRINVLVSTNST HSNLTLKLGVPDSIPHGLISLSSVSQPPAIAYSECFGEDLHGTFKH TVCNRRTDYTLGRI GEIQCPTKADALAVSKRCISSDSIIFSKVHKDSVDCQSSIIDPMT IRNRNKLPSTVGSVT FWPTETSVEAAIPDLASATMLIRLDGYTIQFRSDSNKCSPRFL SLSGCYNCEAGAKLELE HVTDFGTALGILECPSLGYTTYEYVKNLTLEKSIRTMHLNG SHVEAKCYFRCPNSESQTLI RGELIYLFNDDIRHHNQTLSPGLSPKSGSGWDPFGW FKASWLRAIWAILGGTVSLIGVV IYMVFTLCLKVKKS</p>
Purity	Greater than 90% as determined by SDS-PAGE.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C

Tel: 1-631-559-9269 1-516-512-3133

Email: info@creative-biomart.com Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA

for up to one week.

Storage

Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.

Storage Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

GENE INFORMATION

Gene Name

GP

Synonyms

GP; Envelopment polyprotein; M polyprotein

UniProt ID

[P09613](#)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA