

Recombinant Full Length Horse Complement Component C9(C9) Protein, His-Tagged

Cat. No. RFL28089EF Lot. No. (See product label)

SPECIFICATION

Product Overview Recombinant Full Length Horse Complement component C9(C9) Protein (P48770) (22-547aa), fused to N-terminal His tag, was expressed in E. coli.

Species Horse

Source E.coli

ProteinLength Full Length of Mature Protein (22-547)

Form Lyophilized powder

AA Sequence

GPTPNYAPEPEQQSGTPLPIDCRMSSWSEWSECDPCLRQMFRSRSIEVFGQFNGQ
 RCVDA VGDRRQCVPTACEEVEDDCGNDFQCCTGRCIKKRLCNGDNDCGDFSD
 EDDCENDPRPP CRERVVEESELARTAGYGINILGMDPLSTPFDNEYYNGLCDRVRD
 GNTLTYRKPWNLAS LAYETKADKNFRIEHYEQQIQAFRSVIEERRSHFNADFTLKF
 TPTEAKKCKQPEESCNG TDSSENRIFRFAYSKNETYQLFLSYSSKKEKMFLHVKG
 VIQLGKFMRSRDVVLTTTFLD DIKALPTAYEKGEYIAFLETYGTHYSSSGSLGGLYE
 LIYVLDKASMDQKGVELRDIQRCL GFNLDLSLKDKYEVTAKIDKNDCLKRNEKEIVNI
 MDGSLIDDVISLIRGGTRKYAFELKE KLLKGAKTVNVTDFVNWASSLNDAPVLISQRL
 SPIYNLIPVKMKDAHQKKQNLERAIEDY INEFSVRKCHPCQNGGTVIQIDGQCCLCSC
 PIAFEGIACETGKKKIS

Purity Greater than 90% as determined by SDS-PAGE.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Applications	SDS-PAGE
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for multiple use. Avoid repeated freeze-thaw cycles.
Storage Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
GENE INFORMATION	
Gene Name	C9
Synonyms	C9; Complement component C9
UniProt ID	P48770

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA