

Recombinant Full Length *Sodalis Glossinidius* Glycerol-3-Phosphate Acyltransferase(plsY) Protein, His-Tagged

Cat. No. RFL29073SF **Lot. No.** (See product label)

SPECIFICATION

| | |
|-------------------------|---|
| Product Overview | Recombinant Full Length <i>Sodalis glossinidius</i> Glycerol-3-phosphate acyltransferase(plsY) Protein (Q2NWE5) (1-209aa), fused to N-terminal His tag, was expressed in <i>E. coli</i> . |
| Species | <i>Sodalis glossinidius</i> |
| Source | <i>E. coli</i> |
| Protein Length | Full Length (1-209) |
| Form | Lyophilized powder |
| AA Sequence | MNAIAIGMIIFAYLCGSVSNAILICRLARLPDPRTSGSGNPGATNVLRLGGKLAAGV MV FDLVKGMIPVWIGYGLGLPPFWLGLVAIAACLGHYPIFFHFRGGKGVATALGAIA PIGY DLSGLMIGTWLLTALLSGYSSLGAIVSALIAPFYVWWFKPQFTFPVAMLSCVL LRHHDN IQRLWHGQESRIWHRRQKDRDIDRQQSKH |
| Purity | Greater than 90% as determined by SDS-PAGE. |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Storage | Store at -20°C/-80°C upon receipt, aliquoting is necessary for multiple use. Avoid repeated freeze-thaw cycles. |

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Storage Buffer Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Reconstitution We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

GENE INFORMATION

Gene Name plsY

Synonyms plsY; SG0255; Glycerol-3-phosphate acyltransferase; Acyl-PO4 G3P acyltransferase; Acyl-phosphate--glycerol-3-phosphate acyltransferase; G3P acyltransferase; GPAT; Lysophosphatidic acid synthase; LPA synthase

UniProt ID [Q2NWE5](#)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA