

# Recombinant Full Length Prolipoprotein Diacylglyceryl Transferase(Lgt) Protein, His-Tagged

**Cat. No.** RFL29384BF    **Lot. No.** (See product label)

## SPECIFICATION

<b>Product Overview</b>	Recombinant Full Length Prolipoprotein diacylglyceryl transferase(Lgt) Protein (Q7VRF2) (1-284aa), fused to N-terminal His tag, was expressed in E. coli.
<b>Species</b>	Blochmannia floridanus
<b>Source</b>	E.coli
<b>ProteinLength</b>	Full Length (1-284)
<b>Form</b>	Lyophilized powder
<b>AA Sequence</b>	MQFFYTPNPNPIILKVGCSIFYWYGMMYVLSFIFAMWFLKYRKRYIVNDIFSSTSEIEN L LYLNFLGVVIGGRIGYVLFYQWSMLSQDKLWFFKLWEGGMSFHGGLIGVIISIMWF SYCK NQSFLKISDFIVPAVPVGLGLGRLGNFINGELWGRVTFDIPWAMLKFNALLQD LLTLKIH PEWKFFFDYYGALPRHPSQLYEMILEGIVLFVVIYIFSCKSRPEGSISGLFLL LYGLFRI IIEFFRQPDIIHGLINNFITLGQVLSFPMVIFGFIIMYVSYKFK
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Storage</b>	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

<b>Storage Buffer</b>	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
-----------------------	---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**GENE INFORMATION****Gene Name**

lgt

**Synonyms**

lgt; Bfl265; Phosphatidylglycerol--prolipoprotein diacylglyceryl transferase

**UniProt ID**[Q7VRF2](#)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA