

Recombinant COVID-19 Spike S1+S2 (T19R, V70F, E156G, 157-158 deletion, A222V, K417N, L452R, T478K, D614G, P681R, D950N) protein(Met1-Pro1213), His-tagged

Cat. No. Spike-1290V Lot. No. (See product label)

SPECIFICATION

Product Overview	Recombinant SARS-CoV-2 Spike S1+S2 (YP_009724390.1, with mutations T19R, V70F, E156G, 157-158 deletion, A222V, K417N, L452R, T478K, D614G, P681R, D950N) (Met1-Pro1213) was expressed in Insect Cells with a polyhistidine tag at the C-terminus. The mutations were identified in the SARS-CoV-2 variant (known as variant AY.2) which emerged in the India.
Species	Sars-CoV-2
Source	Insect Cells
ProteinLength	Met1-Pro1213
Form	Lyophilized from sterile 20mM Tris, 300mM NaCl, 10% glycerol, pH 8.0. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization.
Molecular Mass	The recombinant SARS-CoV-2 Spike S1+S2 consists of 1207 amino acids and predicts a molecular mass of 134.20 kDa.
Endotoxin	< 1.0 EU per µg protein as determined by the LAL method.
Purity	> 90 % as determined by SDS-PAGE.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Storage Samples are stable for up to twelve months from date of receipt at -20°C to -80°C. Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution It is recommended that sterile water be added to the vial to prepare a stock solution of 0.2 ug/ul. Centrifuge the vial at 4°C before opening to recover the entire contents.

GENE INFORMATION

Gene Name [S surface glycoprotein \[Severe acute respiratory syndrome coronavirus 2 \]](#)

Official Symbol [Spike](#)

Synonyms Spike; S1 protein; Spike glycoprotein Subunit1; S glycoprotein Subunit1; Spike protein S1

Gene ID [43740568](#)

Protein Refseq [YP_009724390.1](#)

UniProt ID [P0DTC2](#)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA