

Dog Uterus Membrane Lysate, Total Protein

Cat. No. Uterus-793D **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Tissues are freshly dissected and homogenized in modified protein lysis buffer supplemented with a cocktail of 7 protease inhibitors to minimize proteolysis.
Species	Dog
Source	Dog Uterus Membrane
Applications	Membrane protein is ready to immediate use for studying membrane proteins such as receptors. Membrane protein is ideal for many downstream applications such as Western blotting, protein assays, enzyme-activity assays, gel electrophoresis, immunoprecipitation, protein-protein interaction analysis, and tissue specific expression identification.
Stability	Protein is stable for one year if stored properly at -70 centigrade.
Storage	Store at -70 centigrade
Concentration	1-2mg/ml

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA