

Yeast Nuclei Isolation Kit

Cat. No. YNI-016K **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Yeast Nuclei Isolation Kit is suitable for extracting complete viable nuclei from various yeasts. The extraction process is simple and convenient. The prepared yeast Nuclei has high purity, maintains natural activity, and has minimal cross-contamination. The extracted nuclei can be used for downstream applications such as nuclear function studies or protein extraction.
Size	50T
Storage	Store at -20°C. Extract B should avoid repeated freezing and thawing. Store the kit as required in the unopened state. The shelf life is one year. Please use the reagent as soon as possible after opening it.
Shipping	Ice pack
Kit Components	Yeast nuclear protein extract A: 20mL, 2-8°C; Yeast nuclear protein extract B: 25mL, -20°C; Yeast nuclear protein extract C: 25mL, 2-8°C; Preservation solution D: 20mL, 2-8°C;
Materials Required but Not Supplied	1. Centrifuge 2. Oscillator 3. Vortex mixer 4. Pipette 5. Refrigerator

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

6. Ice box
7. PBS buffer
8. Centrifuge tube
9. Suction head
10. Disposable gloves

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA