

L-Lactate Colorimetric Assay Kit

Cat. No. Kit-2412 Lot. No. (See product label)

SPECIFICATION

Product Overview

Description

Lactic acid, a major stereoisomer of lactate formed in human intermediary metabolism, has often been used as a physiological indicator for stress. Its determination in serum (normal range in blood 0.5-2.0 mM) is very important in the diagnosis and medical management of various diseases, such as tissue hypoxia, diabetes, circulatory failure and hematological disorders. Unlike a typical lactic acid-LDH assay, this new lactic acid assay kit provides a simple, reliable method for quantifying lactic acid in biological samples such as blood (serum and plasma), culture and fermentation media, etc. In the assay, lactic acid is oxidized by enzyme reactions to yield color product, which can be measured at 570 nm for colorimetric assay, And the color intensity is proportional to lactic acid concentrations, therefore the sample lactic acid concentration can be accurately calculated based on the lactic acid standards.

Applications

For biological research: L-Lactate measurement in biological samples
For drug/pharm research: Drug influence on L-Lactate metabolism

Storage at -80°C

Shipping Icepacks

Size 100 Assays

Kit Components 1. L-Lactate Assay Buffer(10x), 10ml

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

	2. L-Lactate Standard (10x), 100ul 3. Assay Probe, 45ul 4. L-Lactate Enzyme Mix, 55ul
Materials Required but Not Supplied	A plate reader capable of measuring absorbance at 570nm Adjustable pipettes and a repeat pipettor Distilled water(milliQ or HPLC-grade) Clear Flat bottom 96-well plates if not included in the kit purchased
Detection method	Colorimetric method at 570nm
Compatible Sample Types	Plasma,Serum,Cell culture supernatant, other body fluid
Preparation	Reagent Preparation Note: All reagents are frozen. We recommend you spin small vials before opening. 1. L-Lactate Assay Buffer(10x) Note: Please don't do dilution in the provided L-Lactate Assay Buffer bottle itself. The bottle contains 10ml of 10x L-Lactate Assay Buffer. The assay buffer must be equilibrated to room temperature before use. Please make appropriate amount of 1xL-Lactate Assay Buffer for the assay, for example, by diluting 1ml of 10x L-Lactate Assay Buffer with 9 ml of distilled water. Prepare additional diluted assay buffer if needed. Store at -80°C. 2. L-Lactate Standard (10x) Note: Please don't do dilution in the provided L-Lactate standard vial itself. The vial contains 100μl of 1mM L-Lactate Standard. The standard must be equilibrated to room temperature before use. Dilute 100μl of 1mM L-Lactate Standard with 900μl of diluted L-Lactate Assay Buffer to prepare a 100μM L-Lactate Standard. 1ml of diluted standard is enough for making 3 standard curves if assayed in duplicate. Store at -

80°C.

3. Assay Probe

The vial contains 45µl of Assay Probe. It is sufficient for 100 wells. The assay probe must be equilibrated to room temperature before use. Protect from light. Store at -80°C.

4. L-Lactate Enzyme Mix

The vial contains 55µl of L-Lactate Enzyme Mix. The enzyme mix must be thawed on ice before use. It is sufficient for 100 wells. Store at -80°C.

Assay Protocol

Protocol

1. Sample Preparation

Blood, Serum, Plasma, other body fluid, or cell culture supernatant can be measured directly by a series of dilutions of the sample (1/2; 1/4; 1/8; ...) to ensure the readings are within the standard curve range. Your samples can be diluted with dH₂O or diluted L-Lactate Assay Buffer. However, keep in mind if samples (such as hemolyzed serum/plasma) contain high level of lactate dehydrogenase capable of converting lactate to pyruvate, it is important for the samples to be deproteinated.

Add 50µl of samples to each well. We suggest that samples be assayed in duplicate.

Note:

If prepared samples are not assayed the same day, store the samples at -80°C. If samples need to be deproteinated, make sure to deproteinize the samples prior to storing in the freezer. The deproteinated samples will be stable for one month while stored at -80°C. For frozen samples, dilutions of samples must be done right before assaying.

2. Standard Curve Preparation

Add 50 μ l 40 μ l 30 μ l 20 μ l 10 μ l 5 μ l 1 μ l and 0 μ l of diluted L-Lactate Standard to each well, then adjust volume to 50 μ l/well with diluted L-Lactate Assay Buffer.

We suggest that L-Lactate Standards be assayed in duplicate. A standard curve has to be included with samples in each assay.

3. Reaction Solution Preparation

50 μ l Reaction Solution is required for each well, and 1ml of Reaction Solution is prepared with the followings:

982 μ l L-Lactate Assay Buffer(1x)

10 μ l L-Lactate Assay Enzyme Mix

8 μ l Assay Probe

1ml of Reaction Solution is enough for ~20 wells.

Prepare appropriate amount of Reaction Solution for the number of assays to be performed.

Note: Please do not keep any leftover Reaction Solution as it has to be made fresh right before you run assays. If old Reaction Solution is used, assays will not work.

4. Perform the assay

a. Mix the reaction solution well and add 50 μ l of the reaction solution to each well containing the L-Lactate standards and test samples.

b. Incubate the reaction for 30 minutes in a 37°C incubator.

c. Measure absorbance at 570 nm in a microplate reader.

Assay time

30 minutes

Analysis

a. Average the OD570 nm values of replicate wells of each L-Lactate Standard, test samples, and blank. Then, correct background by subtracting the value derived from the blank from all sample readings.

b. Make a standard curve by plotting OD570 nm values from each L-Lactate Standard as a function of L-Lactate concentration. This can be done with excel spreadsheet.

c. Calculate the value of L-Lactate in samples using the equation obtained from the linear regression of the standard curve.

$$\text{L-Lactate}(\mu\text{M}) = [(\text{Corrected absorbance}) - (\text{y-intercept})] / \text{Slope}$$

Sensitivity	2 μ M-100 μ M
--------------------	-------------------------------------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA