

Glutathione Peroxidase Assay Kit

Cat. No. Kit-2449 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Glutathione peroxidase (GPX, EC 1.11.1.9) represents an enzyme family with peroxidase activity whose main biological role is to protect the organism from oxidative damage. It helps prevent lipid peroxidation of cellular membranes by removing free peroxide in the cell. GPX catalyzes the following reaction with glutathione reductase (GR), Simple, direct and high-throughput assays for GPX activity find wide applications. improved assay directly measures NADPH consumption in the enzyme coupled reactions. The measured decrease in optical density at 340 nm is directly proportional to the enzyme activity in the sample.
Storage	-20°C
Shipping	On Ice
Size	100 tests
Detection method	OD340nm
Compatible Sample Types	Biological
Features & Benefits	Sensitive and accurate. Use 10 µL sample. Linear detection range 40 to 800 U/L GPX activity.
Assay time	20 min

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA



Sensitivity

40 U/L

☎ Tel: 1-631-559-9269 1-516-512-3133

☎ Email: info@creative-biomart.com ☎ Fax: 1-631-938-8127

☎ 45-1 Ramsey Road, Shirley, NY 11967, USA