

Recombinant Canine FOLR1 protein, Fc-tagged, APC labeled

Cat. No. FOLR1-1062CA **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Canine FOLR1 protein, Fc-tagged was Conjugated with NH2-Reactive APC(Allophycocyanin).
Species	Canine
Source	HEK293
Description	NH2-Reactive Allophycocyanin is used for antibody conjugation for immunostaining and cellular proteins for tracing. NH2-Reactive Allophycocyanin has succinimidyl ester groups, and can easily make a covalent bond with amino groups of protein without any activation process.
Form	Liquid in sterile PBS buffer pH 7.4 (may contains tiny amount of reaction buffer and Allophycocyanin).
Endotoxin	<1EU/ug by LAL method.
Stability	Sample is stable for at least 12 months when stored at -20 to -80°C from date of manufacture. For longer storage, add equal volume of glycerol to the sample solution and store at -20 to -80°C.
Storage	Store under sterile conditions at -20 to -80°C immediately after receipt. It is recommended that the protein be aliquoted for optimal storage. Avoid freeze-thaw cycles.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Related Product	FOLR1-1062C
------------------------	-------------

Conjugation	APC
--------------------	-----

GENE INFORMATION

Official Symbol	FOLR1
------------------------	-------

Synonyms	FOLR1; folate receptor 1 (adult); folate receptor alpha; folate binding protein 1; Fbp1;
-----------------	--

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA