

Active Recombinant SARS-CoV-2 RBD Protein, His-tagged

Cat. No. RBD-001S **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Purified recombinant fusion protein that includes the RBD of SARS-CoV-2 Spike protein (S-Protein) and Virion Surface Domain of Membrane Protein (M-protein) with C-terminal HIS-tag was expressed in HEK293.
Species	SARS-CoV-2
Source	HEK293
Bio-activity	Binding of SARS-CoV-2 S-M Fusion Protein to Human ACE2: 2 ug/ml of SARS-CoV-2 Spike and Membrane viral fusion protein (gel insert) was immobilized on an ELISA plate and its ability to bind to Human ACE2 was measured. The EC50 was between 10-20 ng/ml.
Molecular Mass	34 kDa
AA Sequence	MALWIDRMQLLSIALSLALVTNSAKGIYQTSNFRVQPTESIVRFPNITNLCPFGEVFN ATRFASVYAWNRRKRISNCVADYSVLYNSASFSTFKCYGVSPTKLNDLCFTNVYADSF VIRGDEVQRQIAPGQTGKIADYNYKL PDDFTGCVIAWNSNNLDSKVGGNYNYLYRLFR KSNLKPFFERDISTEIYQAGSTPCNGVEGFNCYFPLQSYGFQPTNGVGYQPYRVVVL SFELLHAPATVCGPKKSTNLVKNKCVNFNFNGLTGTGVLTESNKKFLGGGSGGGSM ADSNGTITVEELKKHHHHHH
Purity	> 90% as determined by SDS-PAGE and Coomassie blue staining
Stability	Stable for 12 months from the date of receipt of the product under proper storage and

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

handling conditions. Avoid repeated freeze-thaw cycles.

Storage


Store at -80 centigrade.

Concentration

>50 µg/mL as determined by Bradford protein assay method.

Storage Buffer

PBS, 10% glycerol.

**Binding activity
measured by ELISA**



Immobilized SARS-CoV-2 Spike Membrane Fusion Protein (RBD of S-Protein and Virion Surface Domain of M-Protein with His Tag at 2 µg/mL (100 µL/well) can bind to Human ACE2 Fc Tag, the EC₅₀ of the human ACE2 protein is 10-20 ng/mL.

 Tel: 1-631-559-9269 1-516-512-3133


 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

SDS-PAGE



ELISA plates were coated with anti-HIS as capture antibody



5ng of HIS-tagged SARS-Cov-2 Spike-membrane recombinant fusion protein was then added as substrate. Human anti-S antibody against SARS-Cov-2 Spike protein RBD was applied for detection. HRP assay was used for readout. Coating with anti-DDK antibody was served as negative control.